

ViotTM

Video otoscopy **made** effective

Capture
high-quality images
and video for
documentation and
counseling





- Ergonomic design for easy operation and patient comfort.
- Auto white balancing ready for use as soon as the image appears on the
- One button control for a full assessment of both ears without having to revert to the computer to change ear, capture images or video.
- See the ear drum and the lower part of the ear canal simultaneously with 90° angle view.
- Operate through the Viot™ suite or directly from the Affinity2.o/ Equinox2.o, and Callisto™ Suites.
- Windows®10 compatible.
- Full integration to the Noah and OtoAccesss™ databases.



Audiometry Tympanometry ABR OAE Hearing Aid Fitting Balance

2 - 03/2019

Science made smarter

Interacoustics is more than state-ofthe-art solutions

Our mission is clear. We want to lead the way in audiology and balance by translating complexity into clarity:

- Challenges made into clear solutions
- Knowledge made practical
- Invisible medical conditions made tangible and treatable

Our advanced technology and sophisticated solutions ease the lives of healthcare professionals.

We will continue to set the standard for an entire industry. Not for the sake of science. But for the sake of enabling professionals to provide excellent treatment for their millions of patients across the globe.

Interacoustics.com



Related products



Callisto
Portable audiometry,
REM & HIT



Affinity^{2.0}
The integrated fitting solution



Titan Middle ear analyzer

Product specifications

All technical and hardware specifications concerning all products can be downloaded from our website.



ventas@beracahmedica.mx





Audiometry Tympanometry ABR OAE Hearing Aid Fitting Balance