

MODEL: 350100 - CALL CORD WITH TILT/RELEASE DIN



350100



FEATURES

- One (1) momentary push button switch
- Call Assurance LED
- Sealed design keeps infectious contaminants out of the pendant cavity
- Tilt / Release 9-pin DIN connector allows for ease of disconnecting call cord from the patient station
- May be safely washed or sterilized
- Approved for use in oxygen enriched environments

SPECIFICATIONS

Weight: 0.40 lbs (0.18kg)

Size: Length: 10 Ft (3.05m)

Housing & Finish: High-impact molded plastic in light grey with sealed red switch

Terminations: 9-pin DIN connector

DESCRIPTION

The 350100 – Call Cord with Tilt / Release DIN is used in conjunction with Responder 5 patient stations. With a 9-pin DIN connector, the Call Cord plugs securely into the pillow speaker jack on the patient station to provide a unique call-in priority that can be customized based on a facility's needs. When the button on the call cord is pressed then the clear ring around the button illuminates to indicate a call has been placed. The sealed design keeps infectious

contaminants out of the pendant cavity.

The smooth shape and minimal creases allow for easy cleaning.

Distribuidor

Beracah
M E D I C A

✉ ventas@beracahmedica.mx

☎ 662 15 10 979



Rauland is a Division of AMETEK, Inc.

USA +1 800 752 7725

Asia Pacific +65 64835750

www.rauland.com

Mid East & Africa +20 122 2154016

Europe +30 693 7501168

Latin America +1 847 454 4203

Architect and Engineer (A&E) Specifications
available online at customerconnection.rauland.com
Specifications subject to change without notice.